


# Momo

## Features

- Highly manoeuvrable design allows for easy access through narrow doorways and confined areas
- Loop handlebar offers comfortable finger and thumb throttle operation
- Auto-latching lock up mechanism which allows assembly and disassembly in seconds
- Removable swivel seat with folding back making it easy to store in the boot
- Easy to remove battery pack and off-board charger for fast, convenient charging away from scooter
- Contemporary, new stylish design for more comfort and ease of use

## Dimensions cm/in



Overall Dimension : A x B x C 102 x 54 x 82cm/40 x 21 x 32"

Seat Dimension : D x E x F 40 x 38 x 35cm/16 x 15 x 14"

Seat Height from Ground : G 48-53cm/19-21"

Seat Height from Deck : H 37-42cm/15-17"

Wheel Base : I 76cm/30"

## Weight kg/lb

base	battery/ea	seat	total with battery	capacity
27/59	7/15	7/15	47/104	136/300



## Technical Specs

Max Speed	7.2kph/4.5mph
Range	17km/10.6mi
Turn Radius (3 Wheel)	86cm/34"
Turn Radius (4 Wheel)	120cm/47"
Ground clearance	8cm/3"
Motor	DC 24V, 130W (Rating)/480W (Max)
Controller	PG S-Drive 70A
Battery	12V/22A/H x 2pcs
Charger	1.8A off-board
Gradient	6°
Front Wheel	8"(200x50) PU tire
Rear Wheel	220 x 70 PU tire
Brake	electromagnetic brake

**merits**  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.