


Eco 134/144

Features

- Puncture-free tyres are standard
- Loop handlebars offer comfortable finger or thumb throttle operation
- Fully adjustable contoured seat with simple operation of slide, recline and swivel mechanisms to allow maximum user comfort.
- Infinite adjustment of steering tiller means anyone can easily find a comfortable position
- High mount charger socket means you never have to bend over
- Adjustable "coil&shock absorber" suspension that really works and will give superior comfort in all conditions.
- Removable front basket with carry handle
- Adjustable headlight, indicators front and rear, rear tail lights and reverse warning makes this scooter suitable for all occasions

Dimensions cm/in



Overall Dimension (3 Wheel): 121 x 64 x 115cm/47 x 25 x 45"
A x B x C

Overall Dimension (4 Wheel): 122 x 64 x 115cm/48 x 25 x 45"
A x B x C

Seat Dimension : D x E x F 46 x 40 x 51cm/18 x 16 x 20"

Seat Height from Ground : G 56-64cm/22-25"

Seat Height from Deck : H 40-48cm/16-19"

Wheel Base : I 84cm/33"

Weight kg/lb

base	battery/ea	seat	total with battery	capacity
47/104	17/37	18/40	99/218	140/309



Technical Specs

Max Speed	10kph/6.25mph
Range	20-40km/13-27mi
Turn Radius (3 Wheel)	104cm/41"
Turn Radius (4 Wheel)	137cm/54"
Ground clearance	10cm/4"
Motor	DC 24V, 250W (Rating)/915W (Max)
Controller	Rhino 2 90A
Battery	12V/35-50A/H x 2pcs
Charger	5A off-board
Gradient	12°
Front Wheel	10" foam filled tire
Rear Wheel	10" foam filled tire
Brake	intelligent, regenerative, electromagnetic brakes

merits
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.