



Maximum User Weight and Incline Angle

S940A



			10-12km/h
_____	12° (21%)		162kg
_____	10° (18%)		204kg
_____	7.5° (13%)		204kg
_____	0° (0%)		204kg

S941A



			10-12km/h
_____	12° (21%)		126kg
_____	10° (18%)		159kg
_____	7.5° (13%)		159kg
_____	0° (0%)		159kg

S946



			10km/h
_____	12° (21%)		100kg
_____	10° (18%)		136kg
_____	7.5° (13%)		136kg
_____	0° (0%)		136kg

S836/S846



			10-12km/h
_____	12° (21%)		182kg
_____	10° (18%)		182kg
_____	7.5° (13%)		182kg
_____	0° (0%)		182kg

S830/S840

_____	12° (21%)	108kg	10km/h
_____	10° (18%)	136kg	
_____	7.5° (13%)	136kg	
_____	0° (0%)	136kg	

S840 DX

_____	12° (21%)	108kg	10km/h
_____	10° (18%)	136kg	
_____	7.5° (13%)	159kg	
_____	0° (0%)	159kg	

S745

_____	10° (18%)	112kg	10km/h
_____	7.5° (13%)	118kg	
_____	6° (10%)	136kg	
_____	0° (0%)	136kg	

S344

_____	12° (21%)	180kg	10-12km/h
_____	10° (18%)	188kg	
_____	7.5° (13%)	227kg	
_____	0° (0%)	227kg	