

# Regal 344A




Regal 344A XL

## Features

- Loop handlebars offer comfortable finger or thumb throttle operation
- 4 wheel adjustable suspension system allowing for more rider control and comfort
- Large 14" wheels which provide greater ground clearance
- Larger 4 pole motor and larger controller delivering much more power for climbing steep driveways or hills
- Built-in half speed safety switch which operates automatically when turning to provide extra safety
- New improved, fully adjustable leather-feel seat, allowing maximum rider comfort
- Separate emergency braking system for added security
- With new easy to access 12 volt socket
- Available in wide seat version (Regal 344A XL)

## Dimensions cm/in



Overall Dimension : A x B x C	152 x 71 x 135cm/60 x 28 x 53"
Seat Dimension : D x E x F	51 x 46 x 51cm/20 x 18 x 20"
Seat Dimension (Regal 344A XL): D x E x F	80 x 46 x 51cm/31.5 x 18 x 20"
Seat Height from Ground : G	56-64cm/22-25"
Seat Height from Deck : H	40-48cm/16-19"
Wheel Base : I	104cm/41"

## Weight kg/lb

base	battery/ea	seat	total with battery	capacity
81/ 178	25/55	19/42	150/330	227/500



## Technical Specs

Max Speed	12kph/7.5mph
Range	40-50km/26-32mi
Turn Radius	167cm/65.7"
Ground clearance	10cm/4"
Motor	DC 24V, 650W (Rating)/2200W (Max)
Controller	Rhino 160A
Battery	12V/60-80A/H x 2pcs
Charger	8A off-board
Gradient	10°
Front Wheel	14" pneumatic
Rear Wheel	14" pneumatic
Brake	intelligent, regenerative, electromagnetic brakes

**merits**  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.