



Ezi-Go


AUSTRALIA

www.meritsaustralia.com

Ezi-Go




P321XL

Features

- Dismantle in seconds for easy transportation and storage
- Puncture proof tyres front and rear
- Quick release seat and armrests
- Suitable for indoors and shopping centres
- Joystick mounting for left and right hand operation
- Quick release trigger separates the base sections
- Thickly padded seat gives greater comfort to regular users
- Intelligent, regenerative electro magnetic brakes activate automatically when Joy Stick is released
- Available in high back seat version (P321XL)

Dimensions cm



Ezi-Go dismantled

Overall Dimension : A x B x C 86 x 56 x 84cm/34 x 22 x 33"

Seat Dimension (P321) 46 x 38 x 38cm/18 x 15 x 15"
: D x E x F

Seat Dimension (P321XL) 46 x 40 x 46cm/18 x 16 x 18"
: D x E x F

Seat Height from Ground : G 46-53cm/18-21"

Seat Height from Deck : H 33-43cm/13-17"

Weight kg/lb

base	battery/ea	seat	total with battery	capacity
25/55	4.5/10	11/24	45/99	114/250

Technical Specs.

Max Speed	6.4kph/4mph
Range	14.4km/9mi
Turn Radius	77cm/30.3"
Ground clearance	4.5cm/1.8"
Motor	DC 24V 90W(rating)/300W(Max)
Controller	PG VSI 50A
Battery	12V/15A/H x 2pcs
Charger	1.5A off-board
Gradient	6°
Caster	6" x 2" PU tire
Drive Wheel	200 x 75 PU tire
Brake	intelligent, regenerative, electromagnetic brakes

 **merits**
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.