

Ezi-Go




P321XL

Features

- Dismantle in seconds for easy transportation and storage
- Puncture proof tyres front and rear
- Quick release seat and armrests
- Suitable for indoors and shopping centres
- Joystick mounting for left and right hand operation
- Quick release trigger separates the base sections
- Thickly padded seat gives greater comfort to regular users
- Intelligent, regenerative electro magnetic brakes activate automatically when Joy Stick is released
- Available in high back seat version (P321XL)

Dimensions cm



Ezi-Go dismantled

Overall Dimension : A x B x C 86 x 56 x 84cm/34 x 22 x 33"

Seat Dimension (P321) 46 x 38 x 38cm/18 x 15 x 15"
: D x E x F

Seat Dimension (P321XL) 46 x 40 x 46cm/18 x 16 x 18"
: D x E x F

Seat Height from Ground : G 46-53cm/18-21"

Seat Height from Deck : H 33-43cm/13-17"

Weight kg/lb

base	battery/ea	seat	total with battery	capacity
25/55	4.5/10	11/24	45/99	114/250

Technical Specs.

Max Speed	6.4kph/4mph
Range	14.4km/9mi
Turn Radius	77cm/30.3"
Ground clearance	4.5cm/1.8"
Motor	DC 24V 90W(rating)/300W(Max)
Controller	PG VSI 50A
Battery	12V/15A/H x 2pcs
Charger	1.5A off-board
Gradient	6°
Caster	6" x 2" PU tire
Drive Wheel	200 x 75 PU tire
Brake	intelligent, regenerative, electromagnetic brakes

merits
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.