



Fende S946

**merits**  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)

# Fende S946



## Features

- Stylish new modern design with a modesty panel gives this scooter a smart attractive look.
- Adjustable prestige car like suspension front & rear delivers a smooth ride over varied terrain.
- Comfortable leather feel contoured high back seat with slide, swivel, recline & adjustable armrests
- Unique battery storage design allows for greater floor space & a safe secure lockable storage area for personal items.
- New modern LCD dashboard display offers clearer more easy to read battery gauge, speed, time & temperature.
- LED lighting system delivers brighter illumination with less power consumption.
- Built in speed reduction switch which is activated when turning offering exceptional safety.
- With new easy to access 12 volt socket (for phone or GPS)



LED lighting system.



Lockable storage drawer underneath the seat.




Extra storage & Shopping bag hook.



Digital LCD Dashboard Display.

## Dimensions cm/in



Overall Dimension : A x B x C 120 x 63 x 132cm/47.2 x 24.8 x 52"

Seat Dimension : D x E x F 46 x 41 x 46cm/18 x 16 x 18"

Seat Height from Ground : G 64-69cm/25-27"

Seat Height from Deck : H 41.5-46.5cm/16.3-18.3"

Wheel Base : I 88cm/34.6"

## Weight kg/lb

Base	Battery/ea	Seat	Total with Battery	Capacity
70/154	6.5/14	18.5/40.5	114.5/252	136/300

## Technical Specs

Max Speed	10kph/6.25mph
Range	30km/19mi
Turn Radius	143cm/56.3"
Ground clearance	10cm/3.9"
Motor	DC 24V, 295W (Rating)/1150W (Max)
Controller	Rhino 2 90A
Battery	12V/22A/H x 4pcs
Charger	5A off-board
Gradient	12°
Front Wheel	10"
Rear Wheel	10"
Brake	intelligent, regenerative, electromagnetic brakes

**merits**  
**AUSTRALIA**

[www.meritsaustralia.com](http://www.meritsaustralia.com)

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.