



Maverick 10



Maverick 10

Features

- Mid-wheel drive design delivers outstanding maneuverability in tight, compact spaces and a smooth and stable ride both indoors and outdoors even over rough terrain
- Dual in-line motors enhanced efficiency, torque, range and performance
- Full front and rear suspension design offers a comfortable ride
- Low center of gravity design delivers a smooth and stable ride
- Swing back joystick mount is standard
- Height and width adjustable armrest, semi reclining backrest and height adjustable seat for various needs



Excellent stability on footplate



True six wheel design



Flip-up armrest



Height adjustable armrest




Width adjustable armrest



Dynamic Shark Controller

Dimensions cm / in



Overall Dimension : A x B x C	86 x 61 x 123cm / 34 x 24 x 48in
Seat Dimension : D x E x F	46 x 41 x 46cm / 18 x 16 x 18in
Seat Height from Ground : G	51-56cm / 20-22in
Seat Height from Deck : H	43-53cm / 17-21in

Weight kg / lb

base	battery/ea	seat	total with battery	capacity
41 / 90	12 / 26	19 / 42	84kg / 184lb	136kg / 300lb

Technical Specs.

Max Speed	8kph / 5mph
Range	28km / 18mi.
Turn Radius	50cm/ 20"
Ground clearance	5cm / 2"
Motor	DC24V, 160W(Rating) / 520W(Max)
Controller	Dynamic Shark 40A
Battery	12V / 30-40AH x 2pcs
Charger	5A off-board
Gradient	8°
Caster	6" x 2" PU tire
Drive Wheel	10" (260x85) foam filled tire
Brake	intelligent, regenerative, electromagnetic brakes